

SD3

Ultrasonic Pocket Doppler

Ultrasonic Pocket Doppler



SD3
Ultrasonic Pocket Doppler



SD1

Ultrasonic Pocket Doppler

With excellent sensitivity, complete interchangeability and high durability, SD3 series Ultrasonic Pocket Dopplers are ideal for routine fetal heart rate detection by the clinicians.



User-Friendly Design

- Compact design with easy one-hand operation
- Bright OLED screen with real-time FHR display
- Hi-Fi speaker

Interchangeable Waterproof Probe

- Plug-and-play
- 2 MHz probe with deep penetration for later pregnancy and larger patients
- 3 MHz probe detects fetal heart beat as early as 9 weeks
- Automatic power-off while putting the probe back to the rack (SD3 Plus & SD3 Pro)



Work Mode

- Real-time work mode for spot check
- Record & play with audio playback and numeric display

Power Source

- Standard AA batteries for easy replacements
- Built-in Li-ion rechargeable battery ensures long working time
- Unique charger stand and wall mount design



SD3 Lite



SD3



SD3 Plus



SD3 Pro

Configurations

Model	SD3 Lite	SD3	SD3 Plus	SD3 Pro
FHR Display	x	√	√	√
Charger Stand	x	x	o	o
Record & Play	x	x	x	√
AA Batteries	o	o	x	x
Li-ion Battery	x	x	√	√

√ Standard x Not available o Optional